


Conventional calibration

SEC User Training Course



Overview


OMNISEC Training course – Tutorial 4



Why, when and how!


Hardware schematic

Detector schematic




Steps in conventional calibration

Limitations



Software Exercise 3



Discussion of results



Questions


Why, when and how?

- Relatively low Cost
- Easy to use: Pump + Column + Detector
- Conventional used with one Concentration Detector
 - RI (Almost Universal)
 - UV (Second Most Used – Proteins & UV Active)
- Publications usually state Mw '*relative to*' or indicate the use of multiple standards
 - (i.e. range from ~1000Da – 4MDa)
- Column Retention Volume must be known
 - Remember: Separation by Size, not Molecular Weight







Hardware schematic

Conventional calibration system



Detection - RI detector schematic




Differential refraction between solvent and solution results in different signals at the photodiodes.

SEC/GPC - Separation

GPC (also known as size-exclusion chromatography, SEC) has long been used as a key tool for measuring molecular weight

- GPC separates macromolecules in solution according to **size** in a chromatographic column
- After the column, the separated molecules can be analysed by **one or more detectors**




Conventional calibration outline

Conventional GPC is the most widely used calculation method

Concentration
detector: RI
or UV

Series of
standards –
conventional
SEC/GPC


Relative
Molecular
weight



Remember: the columns separate by size not molecular weight so the calibration is only relative

- To measure an unknown sample, the column retention volume must be calibrated in some way
- Use polymer standards of known molecular weight
- Flow rate must be controlled carefully
- Accurate concentration not necessary

A series of standards – 10 to 12 standards



- Series of standards with a range of molecular weights
- Set baselines and limits around each standard to perform calibration (Software exercise 3)

1st step: run standards and create calibration curve

Narrow standards calibration



- Table of Standards on OmniSec V. 5
 - M_p and Vp: molecular weight and retention volume at the peak
- Software exercise 3 – conventional calibration

Conventional Calibration - Homopolymers :CC-0013.vcm

General Column Calibration Peak Detection

Standard Set: [Dropdown]
Curve Fit Order: 5 [Dropdown] Cubic Spline
Vp of Flow Rate Marker: 0
Mw Low: 10 %
Mw High: 10 %
Reset Table


Mp of Mass Marker: 0
Vp of Mass Marker: 0

Del	Use	File Name	Mp	Vp (Std)	Vp FRM	Vp (Corr)
<input type="checkbox"/>	<input checked="" type="checkbox"/>	(18-1.vdt)	8000000	16.0307	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	(17-1.vdt)	1870000	17.1624	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	(18-1.vdt)	940000	18.1054	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	(17-1.vdt)	400000	19.2231	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (170,000; ...	170000	20.3690	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (94,400; 1...	94400	21.1037	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (170,000; ...	46500	22.1827	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (94,400; 1...	18000	23.5545	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (170,000; ...	5450	25.7302	0.000	0.000
<input type="checkbox"/>	<input checked="" type="checkbox"/>	PS Mix (94,400; 1...	2560	27.2191	0.000	0.000

OK Cancel Help

2nd step: run unknown sample

Calibration line and calculation



- Total permeation
 - defines the point at which everything that was injected has passed through the column
- Exclusion limit
 - defines the maximum size of a molecule that can be separated by a column

How the M_w is calculated in conventional calibration

Molecular weight moments

- Number Average (M_n)


$$\overline{M}_n = \frac{\sum c_i}{\sum c_i / M_i}$$

- Weight average (M_w)

$$\overline{M}_w = \frac{\sum c_i M_i}{\sum c_i}$$

- Z-average (M_z)

$$\overline{M}_z = \frac{\sum c_i M_i^2}{\sum c_i M_i}$$



How the M_w is calculated in conventional calibration

Molecular weight moments

M_n = Total mass of material divided by the total number of molecules.

- Mid point of the distribution in terms of numbers of molecules
- Sensitive to low MW species (more molecules in a given mass)

M_w = Multiplying by the molecules mass.

- Weights each chain length according to its weight fraction.
- Mid point of the distribution in terms of polymer weight
- Biased towards larger molecules in the distribution

M_z = Multiplying by the molecules mass again.

- Heavily weighed towards the largest molecules in the sample.
- Sedimentation properties

- Number Average (M_n)

$$\overline{M}_n = \frac{\sum c_i}{\sum c_i / M_i}$$

- Weight average (M_w)



$$\overline{M}_w = \frac{\sum c_i M_i}{\sum c_i}$$

- Z-average (M_z)

$$\overline{M}_z = \frac{\sum c_i M_i^2}{\sum c_i M_i}$$

How the M_w is calculated in conventional calibration

Molecular weight moments



- Number Average (M_n)

$$\overline{M}_n = \frac{\sum c_i}{\sum c_i / M_i}$$

- Weight average (M_w)

$$\overline{M}_w = \frac{\sum c_i M_i}{\sum c_i}$$

- Z-average (M_z)

$$\overline{M}_z = \frac{\sum c_i M_i^2}{\sum c_i M_i}$$


Dispersity Đ

Classification of molecular weight distribution

Type of material	$\text{Đ} = M_w/M_n$
Monodisperse	$= 1.0$
Narrow distribution	< 1.2
Medium distribution	< 2.0
Broad distribution	> 2.0

Dispersity

$$\text{Đ} = \frac{M_w}{M_n}$$



Narrow distribution

$$\text{Đ} < 1.2$$

Broad distribution


$$\text{Đ} > 2.0$$

Advantages of conventional calibration

- Simple **setup**
 - *Only one detector required – RI or UV*
- Accurately known **concentrations** not critical for technique
 - *Approximate concentrations are good enough*
- As a technique – excellent precision (**repeatability**)
 - *Dependent on column and pump performance*

So... where is the disadvantage with conventional calibration?

Overlay of conventional calibration curves



Sample eluting at ret. vol. of 26 ml:

1) PBD calibration curve

- $\text{Log}(M_w) = 3.4$


2) PS calibration curve

- $\text{Log}(M_w) = 3.6$

Each polymer has its own
size to molecular weight
relationship:
 $V_h \approx [\eta] \cdot M$

Conventional calibration of polymers with the same chemistry


Polystyrene (sample) relative to polystyrene (calibration standards)



Sample	M_w	M_n	M_w/M_n
PS235k	231,974	88,527	2.620

Conventional calibration of polymers with different chemistry

Polyethylene (sample) relative to polystyrene (calibration standards)




Conventional calibration results

Sample	M_w	M_n	M_w/M_n
PE1475	164,733	20,326	8.105

Triple Detection method results

Sample	M_w	M_n	M_w/M_n
PE1475	64,129	16,782	3.821


Effect of molecular shape on GPC retention volume



- Size exclusion columns separate by hydrodynamic size and not by M_w
- Therefore, structural differences will affect results:
 - **Conformation**
 - **Branching**

Conventional calibration of polymers with different chemistry

Branched polyethylene (sample) relative to polystyrene (calibration standards)



Conventional calibration results

Sample	M_w	M_n	M_w/M_n
PE1476	171,454	20,314	8.440

Triple Detection method results

Sample	M_w	M_n	M_w/M_n
PE1476	86,508	19,139	4.520

Limitations of conventional calibration

- Every polymer has its own calibration line, which means that M_w values are only accurate for same polymer types
 - Only relative M_w obtained!
 - How close the true and relative M_w values are depends on how close the analyte and standards are in chemical composition and structure
- Any structural change, such as branching, will also affect the accuracy of this value
 - Gives relative M_w even further away from true M_w value!
 - Remember: compare apples with apples or at least a spherical fruit!
- Does not give structural information

Summary

Conventional calibration

- Simple technique to give whole polymer distribution
 - Need to take care with sample/solvent/column compatibility
 - Comparison of samples is easy
- Calibration is main difficulty
 - Data is therefore only relative
- Chromatography conditions need to be carefully controlled
 - Retention volume can be affected by change in conditions
- No structural information
 - Not useful for branched polymers

Software Exercise 3

Conventional Calibration on OmniSEC v5

Objectives

This exercise will instruct you on how to use the OmniSEC v5 software to:

- Estimate the molecular weight of two unknown sample using conventional calibration.

Learning Outcomes

Following this exercise, you will be able to use the OmniSEC v5 software to:

- Set baselines and limits.
- Process conventional calibration data using the OmniSEC software.
- Understand the factors affecting Conventional



Go to exercise

The Triple Detection Method

- Read all the points
- Follow each steps
- Fill in the tables
- Answer questions