

# The OMNISEC™ Triple Detection method explained

## SEC-LS User Training Course




# Overview

## OMNISEC Training course – Tutorial 2




Why do we calibrate?



Steps in the Triple Detection method



Software Exercise 1



Discussion of results




Questions

# How is the data calculated?

## Analysis methods

How do we get from a SEC/GPC chromatograms to the results?



Results by sample and peak	
Parameter	Inj. 3 PS244893 108073 0...
	Peak 1
RV (mL)	16.97
Mn (g/mol)	107,400
Mw (g/mol)	247,100
Mz (g/mol)	433,700
Mw/Mn	2.3
IVw (dL/g)	0.8184
Rh(η)w (nm)	13.9
Rgw (nm)	16.55

Conventional  
calibration

*Tutorial 4*

Universal  
calibration

*Tutorial 5*


Triple detection  
method

*Tutorial 1*



## Hardware schematic

- **Triple Detection:** Homopolymers – RI, LS, Viscometer
- **Tetra Detection:** Copolymers – RI, UV, LS, Viscometer



# GPC System and Detectors



## OMNISEC

### Premium Performance

RI- UV- RALS/LALS - Viscometer

Absolute molecular weight, molecular size, intrinsic viscosity, branching, conformation, protein aggregation.

- Redesigned RALS/LALS with significant sensitivity improvements
- New viscometer design
- Improved RI, including change of position in detector module, and excellent sensitivity



## TDA max

### Benchmark

UV – RALS/LALS – RI - Viscometer




## SEC-MALS 20

### Modular detector

- Multi angle LS – 20 angles
- Absolute molecular weight of proteins, synthetic and natural polymers, as well as molecular size expressed as the radius of gyration,  $R_g$ .

# Summary Triple Detection method

How do we get from a SEC/GPC chromatograms to the results?



Results by sample and peak	
Parameter	Inj. 3 PS244893 108073 0...
	Peak 1
RV (mL)	16.97
Mn (g/mol)	107,400
Mw (g/mol)	247,100
Mz (g/mol)	433,700
Mw/Mn	2.3
IVw (dL/g)	0.8184
Rh( $\eta$ )w (nm)	13.9
Rgw (nm)	16.55

## Run narrow and broad standards

Set baseline and limits  
Create Calculation method (calibration)  
Check Calculation method

## Run samples

Set baseline and limits  
Obtain results

We will learn now why and how we do these steps.





## Why do we calibrate?

- 1 Calculate detector offset
- 2 Determine peak symmetry and band broadening
- 3 Determine calibration constants



1

## Calculate detector offset



Detectors peaks offset is due to:

- the serial configuration in the OMNISEC
- the band broadening due to tubing


Software algorithms take into account these effects.



2

## Determine peak symmetry and band broadening

The symmetry describes the shape of the peaks



Band-broadening looks into:

- The effect of running the sample through the chromatography system
- 'stretching' effect the sample is subjected to.

# Equations Governing The Detectors

## Concentration Detectors

Differential  
Refractive Index

$$RI_i = \frac{K}{n_0} \cdot \left( \frac{dn}{dc} \right) \cdot C_i$$

UV-Vis

$$A_i = K \cdot \frac{dA}{dc} \cdot C_i$$

## Hydrodynamic Size and Intrinsic Viscosity

Four-capillary  
viscometer  
Intrinsic Viscosity

$$\eta_{sp} = \frac{\eta - \eta_0}{\eta_0} = \frac{4\Delta P}{P_i - 2\Delta P} = [\eta] \cdot C$$

## Absolute Molecular Weight Measurement

Static Light Scattering  
Right Angle (90°) and Low Angle (7°)  
Light Scattering

$$\frac{K \cdot \left( \frac{dn}{dc} \right)^2 \cdot C}{R_\theta} = \frac{1}{M_w P(\theta)} + 2A_2C + 3A_3C^2 + \dots$$

## Equations governing the detectors

$$RI \text{ output (mV)} = K_{RI} \cdot dn/dc \cdot \text{concentration}$$

$$UV \text{ output (mV)} = K_{UV} \cdot dA/dc \cdot \text{concentration}$$

$$Visc. \text{ output (mV)} = K_{Visc.} \cdot IV \cdot \text{concentration}$$

$$LS \text{ output (mV)} = K_{LS} \cdot M_w \cdot (dn/dc)^2 \cdot \text{concentration}$$

Unknown

By running the **narrow standard** of known concentration,  $dn/dc$  and  $M_w$  the detectors constants are calculated

$$K_{RI} = \frac{RI \text{ output (mV)}}{dn/dc \cdot \text{concentration}}$$

$$K_{UV} = \frac{UV \text{ output (mV)}}{dA/dc \cdot \text{concentration}}$$


$$K_{Visc.} = \frac{Visc. \text{ output (mV)}}{IV \cdot \text{concentration}}$$

$$K_{LS} = \frac{LS \text{ output (mV)}}{M_w \cdot (dn/dc)^2 \cdot \text{concentration}}$$

# Triple Detection method

## Step-by-step to obtain the $M_w$ of your samples


Calibration with narrow standard (e.g. PS105k in THF)



$$K_{RI} = \frac{RI \text{ output (mV)}}{dn/dc \cdot \text{concentration}}$$

$$K_{LS} = \frac{LS \text{ output (mV)}}{M_w \cdot (dn/dc)^2 \cdot \text{concentration}}$$

Samples analysis – determine your sample  $M_w$



$$dn/dc = \frac{RI \text{ output (mV)}}{K_{RI} \cdot \text{concentration}}$$

$$M_w = \frac{LS \text{ output (mV)}}{K_{LS} \cdot (dn/dc)^2 \cdot \text{concentration}}$$

# Steps in the Triple Detection method?



Set baselines and limits

Calculation method

Calculation method

Detectors

Calibration type

Band broadening

Analysis type

Calibration standard

MALS Parameters

Calibration standard


Name	Mw (g/mol)	IV (dL/g)	Mw/Mn
Polystyrene, 105 kDa, in THF,...	105453	0.48	1.04

Detectors

Active	Detector	Detector constant	Offset (mL)
<input checked="" type="checkbox"/>	RI	9,111,000	0
<input type="checkbox"/>	UV	1,602,000	0
<input checked="" type="checkbox"/>	RALS	3.513E-8	0.1267
<input checked="" type="checkbox"/>	LALS	3.699E-8	0.1267
<input checked="" type="checkbox"/>	DP	1.046	0.21
<input checked="" type="checkbox"/>	IP	0	0.21

Results by sample and peak

Parameter	Inj. 1 PS 244893 108073...
	Peak 1
RV (mL)	16.68
Mn (g/mol)	98,270
Mw (g/mol)	248,100
Mz (g/mol)	446,700
Mw/Mn	2.524
IVw (dL/g)	0.7945
Rh(η)w (nm)	13.74
Rgw (nm)	11.84
M-H a	0.6549
M-H log K (dL/g)	-3.602
RI peak (mV·mL)	290.1
RALS peak (mV·mL)	347.7
LALS peak (mV·mL)	349.9
DP peak (mV·mL)	119.9
MALS peak (mV·mL)	N/C
Calc. dn/dc	N/C
Recoverv (%)	100.3



Results by sample and peak	
Parameter	Inj. 1 PS245k Mw 245245...
	Peak 1
RV (mL)	17.36
Mn (g/mol)	103,700
Mw (g/mol)	235,500
Mz (g/mol)	394,100
Mw/Mn	2.271
IVw (dL/g)	0.8220
Rh(η) <sub>w</sub> (nm)	13.72
Rg <sub>w</sub> (nm)	N/C
M-H a	0.6992
M-H log K (dL/g)	-3.817


## Triple/Tetra Detection GPC/SEC Results

- Absolute molecular weight
- Molecular weight distribution
- Intrinsic viscosity
- Molecular size
- Mark-Houwink coefficients
- % Polymer

# Molecular size by triple detection

## (Random Coil Polymer)

- Two different measurements of size can be performed using Triple Detection GPC
- Hydrodynamic Radius ( $R_h$ )
  - Viscometer and light scattering data
- Radius of Gyration ( $R_g$ )
  - Comparison of light scattering detector intensities
  - Flory-Fox estimation






# Size measurements - $R_h$

## Hydrodynamic Radius ( $R_h$ )

### Triple Detection SEC/GPC – IV and $M_w$

$R_h$  is the radius of an equivalent **solid sphere** that increases the fluid viscosity by the same amount as the macromolecule.


$$[\eta]M = \frac{10}{3} \pi \cdot N_A \cdot R_h^3$$

IV  $M_w$  Constants

### Triple Detection

- Analyze hydrodynamic size from < 1 nm to the exclusion limit of the SEC column (~200 nm)
- No extrapolation or fitting parameters


### Dynamic Light Scattering (DLS) – Zetasizer products

$R_h$  is the radius of an equivalent **sphere that diffuses** with the same speed as the molecule of interest.

# Size measurements - $R_g$

## Radius of Gyration ( $R_g$ )

$R_g$  is the root-mean-square of the radii from the centre of the mass to the different mass cores within the molecule.



- Direct measurement by changes in scattered light intensities with observation angle
  - RALS/LALS
  - MALS

### Limitations:

- Requires good S/N light scattering signal
- Lower size detection limit = 10-15 nm
  - Limit of Anisotropic scattering
- Large structures require non-linear curve fitting

# Size measurements - $R_g$

## Radius of Gyration ( $R_g$ )




A good estimate of  $R_g$  can be made using viscometry and the **Flory-Fox equation**

$$R_g = \frac{1}{\sqrt{6}} \left( \frac{[\eta] M_w}{\Phi} \right)^{\frac{1}{3}}$$


Equation relates the intrinsic viscosity of a **flexible coil molecule** in solution in terms of  $R_g$

- Where  $\Phi$  = Flory's universal constant
- Only OmniSEC v5

# Comparison of $R_h$ and $R_g$ distribution for broad polystyrene 235k



- Triple detection technique uses 3 detectors:
- Specific equations govern each detector which allow
  - concentration
  - $M_w$
  - intrinsic viscosity
  - molecular sizeto be determined across the entire distribution.
- Adding a fourth detector – UV – allows for compositional analysis in copolymer samples.



## Summary



Questions?



# Fale comigo!

## Gustavo G. Parra, PhD



• **Cel:** +55 11 94556-5914

**E-mail:** [gustavo.parra@pensabio.com.br](mailto:gustavo.parra@pensabio.com.br)